

For Research Use Only

mouse ADCY8-Specific Polyclonal antibody



Catalog Number: 19773-1-AP

Basic Information

Catalog Number:

19773-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009623

GeneID (NCBI):

11514

UNIPROT ID:

P97490

Full Name:

adenylate cyclase 8

Calculated MW:

140 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

ADCY8 belongs to the adenylyl cyclase class-4/guanylyl cyclase family. ADCY8 is a membrane-bound, calcium-stimulable adenylyl cyclase. It may be involved in learning, in memory and in drug dependence. ADCY8 catalyze the reaction: ATP = 3',5'-cyclic AMP + diphosphate. This antibody is specific to ADCY8.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

