

For Research Use Only

ZBTB37 Polyclonal antibody

Catalog Number: 19817-1-AP



Basic Information

Catalog Number:

19817-1-AP

Size:

147 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13864

GenBank Accession Number:

BC003116

GeneID (NCBI):

84614

UNIPROT ID:

Q5TC79

Full Name:

zinc finger and BTB domain
containing 37

Calculated MW:

503 aa, 56 kDa

Observed MW:

56-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 19817-1-AP (ZBTB37 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.