

For Research Use Only

ZNF235 Polyclonal antibody

Catalog Number: 19818-1-AP

1 Publications



Basic Information

Catalog Number:

19818-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13866

GenBank Accession Number:

BC002663

GeneID (NCBI):

9310

UNIPROT ID:

Q14590

Full Name:

zinc finger protein 235

Calculated MW:

738 aa, 84 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF, IHC

Species Specificity:

human

Cited Species:

rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xinliang Zou	35900670	J Cardiovasc Transl Res	IF, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

