

For Research Use Only

GULP1 Polyclonal antibody

Catalog Number: 19902-1-AP **1 Publications**



Basic Information

Catalog Number:

19902-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13829

GenBank Accession Number:

BC001103

GeneID (NCBI):

51454

UNIPROT ID:

Q9UBP9

Full Name:

GULP, engulfment adaptor PTB domain containing 1

Calculated MW:

167 aa, 19 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells, NIH/3T3 cells, HepG2 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Dong Yan	37573356	Cell Death Dis	WB,IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19902-1-AP (GULP1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.