

For Research Use Only

# NUPL1 Polyclonal antibody

Catalog Number: 19907-1-AP

2 Publications



## Basic Information

Catalog Number:

19907-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13901

GenBank Accession Number:

BC001104

GeneID (NCBI):

9818

UNIPROT ID:

Q9BVL2

Full Name:

nucleoporin like 1

Calculated MW:

599 aa, 61 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HepG2 cells, HeLa cells, L02 cells, mouse brain tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Melissa Kane	30084827	Elife	WB
Fani Panagaki	39036650	Nat Phys	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 19907-1-AP (NUPL1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.