

For Research Use Only

# TMEM9 Polyclonal antibody

Catalog Number: 19918-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

19918-1-AP

Size:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13831

GenBank Accession Number:

BC001106

GeneID (NCBI):

252839

UNIPROT ID:

Q9POT7

Full Name:

transmembrane protein 9

Calculated MW:

21 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : LNCaP cells,

## Background Information

The function of TMEM9 remains largely unknown. TMEM9 belongs to the TMEM9 family, and may be involved in intracellular transport.

## Notable Publications

Author	Pubmed ID	Journal	Application
Jianchao He	36596527	Biol Pharm Bull	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



LNCaP cells were subjected to SDS PAGE followed by western blot with 19918-1-AP (TMEM9 antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.