

For Research Use Only

# TMEM9 Polyclonal antibody

Catalog Number: 19918-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 19918-1-AP	<b>GenBank Accession Number:</b> BC001106	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 252839	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9P0T7	
<b>Isotype:</b> IgG	<b>Full Name:</b> transmembrane protein 9	
<b>Immunogen Catalog Number:</b> AG13831	<b>Calculated MW:</b> 21 kDa	
	<b>Observed MW:</b> 23 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : LNCaP cells,
<b>Cited Applications:</b> WB	
<b>Species Specificity:</b> human, mouse, rat	
<b>Cited Species:</b> human	

## Background Information

The function of TMEM9 remains largely unknown. TMEM9 belongs to the TMEM9 family, and may be involved in intracellular transport.

## Notable Publications

Author	Pubmed ID	Journal	Application
Jianchao He	36596527	Biol Pharm Bull	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



LNCaP cells were subjected to SDS PAGE followed by western blot with 19918-1-AP (TMEM9 antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours.