

For Research Use Only

FSHR-Specific Polyclonal antibody

Catalog Number: 19948-1-AP



Basic Information

Catalog Number:

19948-1-AP

Size:

170 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000145

GeneID (NCBI):

2492

UNIPROT ID:

P23945

Full Name:

follicle stimulating hormone receptor

Calculated MW:

78 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

FSHR, also named as LGR1, belongs to the G-protein coupled receptor 1 family and FSH/LSH/TSH subfamily. FSHR is a receptor for follicle-stimulating hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase. Defects in FSHR are a cause of ovarian dysgenesis type 1 (ODG1). Defects in FSHR are a cause of ovarian hyperstimulation syndrome (OHSS). This antibody is specific to FSHR.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

