

For Research Use Only

Huntingtin Polyclonal antibody

Catalog Number: 19957-1-AP



Basic Information

Catalog Number: 19957-1-AP	GenBank Accession Number: NM_002111	Purification Method: Antigen affinity purification
Size: 267 µg/ml	GeneID (NCBI): 3064	Recommended Dilutions: IHC 1:50-1:500 IF 1:50-1:500
Source: Rabbit	UNIPROT ID: P42858	
Isotype: IgG	Full Name: huntingtin	
	Calculated MW: 348 kDa	

Applications

Tested Applications: IF/ICC, IHC, ELISA	Positive Controls: IHC : mouse brain tissue, IF : SH-SY5Y cells,
Species Specificity: human, mouse	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

HTT, also named as HD and IT15, belongs to the huntingtin family. HTT may play a role in microtubule-mediated transport or vesicle function. Defects in HTT are the cause of Huntington disease (HD) which is an autosomal dominant neurodegenerative disorder characterized by involuntary movements (chorea), general motor impairment, psychiatric disorders and dementia. The antibody recognizes the C-term of HTT.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 19957-1-AP (Huntingtin Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 19957-1-AP (Huntingtin Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed SH-SY5Y cells using 19957-1-AP (Huntingtin antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).