

For Research Use Only

HBD-Specific Polyclonal antibody

Catalog Number: 19992-1-AP



Basic Information

Catalog Number:

19992-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000519

GeneID (NCBI):

3045

UNIPROT ID:

P02042

Full Name:

hemoglobin, delta

Calculated MW:

16 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

HBD, also named as Delta-globin, belongs to the globin family. It is involved in oxygen transport from the lung to the various peripheral tissues. This antibody is specific to HBD. It has no cross reaction with HBB.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

