

For Research Use Only

HBQ1-Specific Polyclonal antibody

Catalog Number: 19997-1-AP

1 Publications



Basic Information

Catalog Number:

19997-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005331

GeneID (NCBI):

3049

UNIPROT ID:

P09105

Full Name:

hemoglobin, theta 1

Calculated MW:

16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : recombinant protein,

Background Information

HBQ1, also named as Theta-1-globin, belongs to the globin family. The antibody is specific to HBQ1.

Notable Publications

Author	Pubmed ID	Journal	Application
Kyungho Kim	38067210	Cancers (Basel)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB results of 19997-1-AP (HBQ1 antibody) at dilution of 1:2000.