

For Research Use Only

FZD5 Polyclonal antibody

Catalog Number: 20007-1-AP



Basic Information

Catalog Number:

20007-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003468

GeneID (NCBI):

7855

UNIPROT ID:

Q13467

Full Name:

frizzled homolog 5 (Drosophila)

Calculated MW:

65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

FZD5, also named as HFZ5, belongs to the G-protein coupled receptor Fz/Smo family. It is a receptor for Wnt proteins. FZD5 may be involved in transduction and intercellular transmission of polarity information during tissue morphogenesis and/or in differentiated tissues. FZD5 interacts specifically with Wnt5A to induce the beta-catenin pathway.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data