

For Research Use Only

# AKAP12 Polyclonal antibody

Catalog Number: 20016-1-AP



## Basic Information

Catalog Number:

20016-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005100

GeneID (NCBI):

9590

UNIPROT ID:

Q02952

Full Name:

A kinase (PRKA) anchor protein 12

Calculated MW:

191 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

AKAP12, also named as AKAP250, is an anchoring protein that mediates the subcellular compartmentation of protein kinase A (PKA) and protein kinase C (PKC). The antibody recognizes the C-term of AKAP12.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

