

For Research Use Only

ARHGEF17 Polyclonal antibody

Catalog Number: 20021-1-AP



Basic Information

Catalog Number:

20021-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014786

GeneID (NCBI):

9828

UNIPROT ID:

Q96PE2

Full Name:

Rho guanine nucleotide exchange factor (GEF) 17

Calculated MW:

222 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

ARHGEF17, also named as KIAA0337, TEM4 and p164RhoGEF, acts as guanine nucleotide exchange factor (GEF) for RhoA GTPases. The antibody recognizes the C-term of ARHGEF17.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

