

# AKAP11-Specific Polyclonal antibody

Catalog Number: 20034-1-AP

## Basic Information

**Catalog Number:**

20034-1-AP

**Size:**

210 µg/ml

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_016248

**GeneID (NCBI):**

11215

**UNIPROT ID:**

Q9UKA4

**Full Name:**

A kinase (PRKA) anchor protein 11

**Calculated MW:**

211 kDa

**Observed MW:**

220 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human

## Background Information

AKAP11, also named as AKAP220, KIAA0629 and PRKA11, belongs to the AKAP110 family. AKAP11 binds to type II regulatory subunits of protein kinase A and anchors/targets them. The antibody is specific to AKAP11.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

