

For Research Use Only

# GPR124 Polyclonal antibody, PBS Only

Catalog Number:20050-1-PBS



## Basic Information

Catalog Number:

20050-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_032777

GeneID (NCBI):

25960

UNIPROT ID:

Q96PE1

Full Name:

G protein-coupled receptor 124

Calculated MW:

143 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

GPR124, also named KIAA1531 and TEM5, belongs to the G-protein coupled receptor 2 family and LN-TM7 subfamily. GPR124 is an Orphan receptor. GPR124 is a Tumor endothelial marker (TEM) that is significantly up-regulated during angiogenesis and neoangiogenesis that are crucial for the growth of solid tumors localized on the cell surface and conserved across species are of particular interest for future development of anti-angiogenic therapies. GPR124 is a member of the adhesion family of G protein-coupled receptors and is localized on the surface of endothelial cells. GPR124 is a seven-pass transmembrane receptor, unlike TEM1, TEM7, and TEM8 which span the membrane once. Proteolytically processed soluble GPR124 mediates endothelial cell survival during angiogenesis by linking integrin to glycosaminoglycans. The antibody recognizes near the C-term of GPR124. The antibody can't recognize the isoform3 of GPR124.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse small intestine tissue slide using 20050-1-AP (GPR124 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20050-1-PBS in a different storage buffer formulation.

Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using 20050-1-AP (GPR124 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20050-1-PBS in a different storage buffer formulation.