

NOM1 Polyclonal antibody

Catalog Number: 20066-1-AP

Basic Information

Catalog Number:

20066-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_138400

GeneID (NCBI):

64434

UNIPROT ID:

Q5C9Z4

Full Name:

nucleolar protein with MIF4G domain
1

Calculated MW:

96 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

NOM1, also named as C7orf3, belongs to the CWC22 family. NOM1 plays a role in targeting PPP1CA to the nucleolus (Pubmed:17965019). NOM1 is a new member of the MIF4G/MA3 family of proteins and suggest that it may provide an essential function in metazoans (Pubmed:15715967). The antibody recognizes the C-term of NOM1.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

