

For Research Use Only

NOM1 Polyclonal antibody

Catalog Number: 20066-1-AP



Basic Information

Catalog Number: 20066-1-AP	GenBank Accession Number: NM_138400	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 64434	
Source: Rabbit	UNIPROT ID: Q5C9Z4	
Isotype: IgG	Full Name: nucleolar protein with MIF4G domain 1	
	Calculated MW: 96 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

NOM1, also named as C7orf3, belongs to the CWC22 family. NOM1 plays a role in targeting PPP1CA to the nucleolus (Pubmed:17965019). NOM1 is a new member of the MIF4G/MA3 family of proteins and suggest that it may provide an essential function in metazoans (Pubmed:15715967). The antibody recognizes the C-term of NOM1.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data