

For Research Use Only

SH3TC2 Polyclonal antibody, PBS Only

Catalog Number:20067-1-PBS



Basic Information

Catalog Number:

20067-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_024577

GeneID (NCBI):

79628

UNIPROT ID:

Q8TF17

Full Name:

SH3 domain and tetratricopeptide repeats 2

Calculated MW:

145 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

SH3TC2, also named as KIAA1985, contains 1 SH3 domain and TPR repeats. Defects in SH3TC2 are the cause of Charcot-Marie-Tooth disease type 4C (CMT4C) [MIM:601596]. The antibody can recognize all the isoforms except isoform4 of SH3TC2.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data