

For Research Use Only

# PCDH11Y Polyclonal antibody

Catalog Number: 20069-1-AP



## Basic Information

Catalog Number:

20069-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_032973

GeneID (NCBI):

83259

UNIPROT ID:

Q9BZA8

Full Name:

protocadherin 11 Y-linked

Calculated MW:

147 kDa

Observed MW:

114-147 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, HeLa cells

## Background Information

PCDH11Y, also named as PCDH11, PCDH22 and PCDHY, is a potential calcium-dependent cell-adhesion protein. The antibody can recognize PCDH11X and PCDH11Y.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 20069-1-AP (PCDH11Y antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.