

For Research Use Only

BAGE3 Polyclonal antibody

Catalog Number: 20072-1-AP



Basic Information

Catalog Number:

20072-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_182481

GeneID (NCBI):

85318

UNIPROT ID:

Q86Y29

Full Name:

B melanoma antigen family, member 3

Calculated MW:

12 kDa

Observed MW:

12 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-13 cells,

Background Information

BAGE3, also named as CT2.3, belongs to the BAGE family. BAGE3 is a candidate gene encoding tumor antigens. It is not expressed in normal tissues except in testis, and also expressed in melanoma, bladder and lung carcinomas.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-13 cells were subjected to SDS PAGE followed by western blot with 20072-1-AP (BAGE3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.