

For Research Use Only

SCARF2-Specific Polyclonal antibody

Catalog Number: 20073-1-AP



Basic Information

Catalog Number:

20073-1-AP

Size:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_182895

GeneID (NCBI):

91179

UNIPROT ID:

Q96GP6

Full Name:

scavenger receptor class F, member 2

Calculated MW:

92 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

SCARF2, also named as SREC2 and SREPCR, is probable adhesion protein, which mediates homophilic and heterophilic interactions. In contrast to SCARF1, it poorly mediates the binding and degradation of acetylated low density lipoprotein (Ac-LDL). This antibody is specific to SCARF2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

