

For Research Use Only

ALPK2-Specific Polyclonal antibody

Catalog Number: 20079-1-AP



Basic Information

Catalog Number:

20079-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_052947

GeneID (NCBI):

115701

UNIPROT ID:

Q86TB3

Full Name:

alpha-kinase 2

Calculated MW:

237 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

ALPK2, also named as HAK, belongs to the protein kinase superfamily, Alpha-type protein kinase family and ALPK subfamily. ALPK2 is a kinase that recognizes phosphorylation sites in which the surrounding peptides have an alpha-helical conformation. This antibody is specific to ALPK2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

