

For Research Use Only

CDYL2-Specific Polyclonal antibody

Catalog Number: 20084-1-AP



Basic Information

Catalog Number:

20084-1-AP

Size:

360 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_152342

GeneID (NCBI):

124359

UNIPROT ID:

Q8N8U2

Full Name:

chromodomain protein, Y-like 2

Calculated MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

CDYL2 acts as repressor of transcription. The antibody is specific to CDYL2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

