

For Research Use Only

# CSTL1-Specific Polyclonal antibody, PBS Only

Catalog Number: 20086-1-PBS



## Basic Information

Catalog Number:

20086-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_138283

GeneID (NCBI):

128817

UNIPROT ID:

Q9H114

Full Name:

cystatin-like 1

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse

## Background Information

CSTL1, also named as RCET11, belongs to the cystatin family. The antibody is specific to CSTL1.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data