

For Research Use Only

LDLRAD2-Specific Polyclonal antibody



Catalog Number: 20095-1-AP

Basic Information

Catalog Number: 20095-1-AP	GenBank Accession Number: NM_001013693	Purification Method: Antigen affinity purification
Size: 280 µg/ml	GeneID (NCBI): 401944	
Source: Rabbit	UNIPROT ID: Q5SZ11	
Isotype: IgG	Full Name: low density lipoprotein receptor class A domain containing 2	
	Calculated MW: 29 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

LDLRAD2 belongs to the LDLR family. This antibody is specific to LDLRAD2.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data