

For Research Use Only

LDLRAD2-Specific Polyclonal antibody



Catalog Number: 20095-1-AP

Basic Information

Catalog Number:

20095-1-AP

Size:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001013693

GeneID (NCBI):

401944

UNIPROT ID:

Q5SZ11

Full Name:

low density lipoprotein receptor class

A domain containing 2

Calculated MW:

29 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

LDLRAD2 belongs to the LDLR family. This antibody is specific to LDLRAD2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

