

For Research Use Only

LDLRAD1 Polyclonal antibody

Catalog Number: 20096-1-AP



Basic Information

Catalog Number:

20096-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001010978

GeneID (NCBI):

388633

UNIPROT ID:

Q5T700

Full Name:

low density lipoprotein receptor class
A domain containing 1

Calculated MW:

22 kDa

Observed MW:

13 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue, mouse liver tissue

Background Information

LDLRAD1 Belongs to the LDLR family. The antibody can recognize all the 3 isoforms of LDLRAD1.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 20096-1-AP (LDLRAD1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.