

For Research Use Only

NBPF4-Specific Polyclonal antibody

Catalog Number:20099-1-AP



Basic Information

Catalog Number:

20099-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001143989

GeneID (NCBI):

148545

UNIPROT ID:

Q96M43

Full Name:

neuroblastoma breakpoint family,
member 4

Calculated MW:

72 kDa

Observed MW:

80-90 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

NBPF4 belongs to the NBPF family. This antibody is specific to NBPF4.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20099-1-AP (NBPF4-Specific antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.