

For Research Use Only

# NBPF4-Specific Polyclonal antibody, PBS Only

Catalog Number:20099-1-PBS



## Basic Information

Catalog Number:

20099-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001143989

GeneID (NCBI):

148545

UNIPROT ID:

Q96M43

Full Name:

neuroblastoma breakpoint family,  
member 4

Calculated MW:

72 kDa

Observed MW:

80-90 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

NBPF4 belongs to the NBPF family. This antibody is specific to NBPF4.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20099-1-AP (NBPF 4-Specific antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20099-1-PBS in a different storage buffer formulation.