

ZNF600-Specific Polyclonal antibody

Catalog Number: 20100-1-AP

Basic Information

Catalog Number:

20100-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_198457

GeneID (NCBI):

162966

UNIPROT ID:

Q6ZNG1

Full Name:

zinc finger protein 600

Calculated MW:

83 kDa

Observed MW:

83 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue,

Background Information

ZNF600 belongs to the krueppel C2H2-type zinc-finger protein family. ZNF600 may be involved in transcriptional regulation. This antibody is specific to ZNF600.

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 20100-1-AP (ZNF600-Specific Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.