

For Research Use Only

ARHGEF40 Polyclonal antibody, PBS Only

Catalog Number: 20110-1-PBS



Basic Information

Catalog Number: 20110-1-PBS	GenBank Accession Number: BC000084	Purification Method: Antigen affinity purification
Size: 1 mg/ml	GeneID (NCBI): 55701	
Source: Rabbit	UNIPROT ID: Q8TER5	
Isotype: IgG	Full Name: hypothetical protein FLJ10357	
Immunogen Catalog Number: AG13886	Calculated MW: 414 aa, 46 kDa	
	Observed MW: 150-160 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data