

For Research Use Only

WDR13 Polyclonal antibody

Catalog Number: 20122-1-AP

1 Publications



Basic Information

Catalog Number:

20122-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13793

GenBank Accession Number:

BC002507

GeneID (NCBI):

64743

UNIPROT ID:

Q9H1Z4

Full Name:

WD repeat domain 13

Calculated MW:

485 aa, 54 kDa

Observed MW:

43 kDa, 53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

bovine

Positive Controls:

WB : mouse lung tissue, human heart tissue, Jurkat cells, mouse kidney tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yuying Fu	31050064	Cell Biol Int	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 20122-1-AP (WDR13 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.