

For Research Use Only

# WDR13 Polyclonal antibody

Catalog Number: 20122-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 20122-1-AP	<b>GenBank Accession Number:</b> BC002507	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 300 µg/ml	<b>GeneID (NCBI):</b> 64743	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9H1Z4	
<b>Isotype:</b> IgG	<b>Full Name:</b> WD repeat domain 13	
<b>Immunogen Catalog Number:</b> AG13793	<b>Calculated MW:</b> 485 aa, 54 kDa	
	<b>Observed MW:</b> 43 kDa, 53 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse lung tissue, human heart tissue, Jurkat cells, mouse kidney tissue
<b>Cited Applications:</b> WB, IF	
<b>Species Specificity:</b> human, mouse, rat	
<b>Cited Species:</b> bovine	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yuying Fu	31050064	Cell Biol Int	WB,IF

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 20122-1-AP (WDR13 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.