

TMEM214 Polyclonal antibody

Catalog Number: 20125-1-AP

1 Publications

Basic Information

Catalog Number:

20125-1-AP

Size:

167 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13884

GenBank Accession Number:

BC002467

GeneID (NCBI):

54867

UNIPROT ID:

Q6NUQ4

Full Name:

transmembrane protein 214

Calculated MW:

689 aa, 77 kDa

Observed MW:

77 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IF 1:25-1:100

Applications

Tested Applications:

IF, WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : C2C12,

IF : human muscle slides,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Phu Le Thanh	28214269	Neuromuscul Disord	WB,IF

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



WB result of TMEM214 antibody (20125-1-AP) with C2C12 cells (Control and KD).



Human Muscle Section

IF result of TMEM214 antibody (20125-1-AP) with human muscle slides.