

For Research Use Only

TMEM214 Polyclonal antibody

Catalog Number: 20125-1-AP **1 Publications**



Basic Information

Catalog Number: 20125-1-AP	GenBank Accession Number: BC002467	Purification Method: Antigen affinity purification
Size: 167 µg/ml	GeneID (NCBI): 54867	Recommended Dilutions: WB 1:500-1:1000 IF 1:25-1:100
Source: Rabbit	UNIPROT ID: Q6NUQ4	
Isotype: IgG	Full Name: transmembrane protein 214	
Immunogen Catalog Number: AG13884	Calculated MW: 689 aa, 77 kDa	
	Observed MW: 77 kDa	

Applications

Tested Applications: IF, WB, ELISA	Positive Controls: WB : C2C12, IF : human muscle slides,
Cited Applications: WB, IF	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Phu Le Thanh	28214269	Neuromuscul Disord	WB,IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of TMEM214 antibody (20125-1-AP) with C2C12 cells (Control and KD).



Human Muscle Section

IF result of TMEM214 antibody (20125-1-AP) with human muscle slides.