

For Research Use Only

FMNL1 Polyclonal antibody

Catalog Number: 20152-1-AP **1 Publications**



Basic Information

Catalog Number: 20152-1-AP	GenBank Accession Number: BC021906	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 752	
Source: Rabbit	UNIPROT ID: O95466	
Isotype: IgG	Full Name: formin-like 1	
Immunogen Catalog Number: AG13697	Calculated MW: 122 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB, IHC

Species Specificity:
human, mouse, rat

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Gui Ma	37454877	Mol Cell Probes	IHC,WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data