

PLCXD1 Polyclonal antibody

Catalog Number: 20161-1-AP

Basic Information

Catalog Number:

20161-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14086

GenBank Accession Number:

BC005028

GeneID (NCBI):

55344

UNIPROT ID:

Q9NUJ7

Full Name:

phosphatidylinositol-specific
phospholipase C, X domain
containing 1

Calculated MW:

323 aa, 37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

