

# TMEM121 Polyclonal antibody

Catalog Number: 20165-1-AP

## Basic Information

**Catalog Number:**

20165-1-AP

**Size:**

N/A

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG13891

**GenBank Accession Number:**

BC004221

**GeneID (NCBI):**

80757

**UNIPROT ID:**

Q9BTD3

**Full Name:**

transmembrane protein 121

**Calculated MW:**

319 aa, 36 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

human, mouse, rat

## Background Information

### Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

