

For Research Use Only

TBC1D9B Polyclonal antibody

Catalog Number: 20177-1-AP



Basic Information

Catalog Number: 20177-1-AP	GenBank Accession Number: BC000701	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 23061	
Source: Rabbit	UNIPROT ID: Q66K14	
Isotype: IgG	Full Name: TBC1 domain family, member 9B (with GRAM domain)	
Immunogen Catalog Number: AG14019	Calculated MW: 1250 aa, 140 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data