

For Research Use Only

C20orf4 Polyclonal antibody

Catalog Number: 20179-1-AP



Basic Information

Catalog Number:

20179-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14039

GenBank Accession Number:

BC001751

GeneID (NCBI):

25980

UNIPROT ID:

Q9Y312

Full Name:

chromosome 20 open reading frame 4

Calculated MW:

384 aa, 44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : MCF-7 cells, HeLa cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 20179-1-AP (C20orf4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.