

For Research Use Only

# C20orf4 Polyclonal antibody

Catalog Number: 20179-1-AP



## Basic Information

<b>Catalog Number:</b> 20179-1-AP	<b>GenBank Accession Number:</b> BC001751	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 600 µg/ml	<b>GeneID (NCBI):</b> 25980	<b>Recommended Dilutions:</b> WB 1:200-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9Y312	
<b>Isotype:</b> IgG	<b>Full Name:</b> chromosome 20 open reading frame 4	
<b>Immunogen Catalog Number:</b> AG14039	<b>Calculated MW:</b> 384 aa, 44 kDa	
	<b>Observed MW:</b> 44 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : MCF-7 cells, HeLa cells, HepG2 cells
<b>Species Specificity:</b> human, mouse, rat	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 20179-1-AP (C20orf4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.