

For Research Use Only

POTEG-Specific Polyclonal antibody

Catalog Number: 20202-1-AP



Basic Information

Catalog Number:

20202-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001005356

GeneID (NCBI):

404785

UNIPROT ID:

Q6S5H5

Full Name:

POTE ankyrin domain family,
member G

Calculated MW:

57 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

POTEG, also named as A26C2, POTE14, belongs to the POTE family. POTE has an important signaling function in the reproductive system. This antibody is specific to POTEG.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

