

For Research Use Only

# POTEB-Specific Polyclonal antibody

Catalog Number: 20205-1-AP



## Basic Information

Catalog Number:

20205-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_207355

GeneID (NCBI):

339010

UNIPROT ID:

Q6S5H4

Full Name:

POTE ankyrin domain family,  
member B

Calculated MW:

66 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

POTEB, also named as A26B1 and POTE15, belongs to the POTE family. POTE has an important signaling function in the reproductive system. This antibody is specific to POTEB.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

