

For Research Use Only

# NPAS4-Specific Polyclonal antibody

Catalog Number: 20208-1-AP



## Basic Information

Catalog Number:

20208-1-AP

Size:

130  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_178864

GeneID (NCBI):

266743

UNIPROT ID:

Q8IUM7

Full Name:

neuronal PAS domain protein 4

Calculated MW:

87 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

NPAS4, also named as NXF, acts as a transcriptional activator in the presence of ARNT. NPAS4 can activate the CME (CNS midline enhancer) element and the expression of the drebrin gene. [PMID:16956324] This antibody is specific to NPAS4.

## Storage

Storage:

Store at  $-20^{\circ}\text{C}$ . Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for  $-20^{\circ}\text{C}$  storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

