

For Research Use Only

GPRIN3-Specific Polyclonal antibody, PBS Only

Catalog Number: 20212-1-PBS



Basic Information

Catalog Number:

20212-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_198281

GeneID (NCBI):

285513

UNIPROT ID:

Q6ZVF9

Full Name:

GPRIN family member 3

Calculated MW:

82 kDa

Observed MW:

82 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

GPRIN3, also named as KIAA2027, may be involved in neurite outgrowth. The antibody is specific to GPRIN3.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 20212-1-AP (GPRIN3-Specific antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20212-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 20212-1-AP (GPRIN3-Specific antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 20212-1-PBS in a different storage buffer formulation.