

For Research Use Only

SLIT3-Specific Polyclonal antibody

Catalog Number: 20218-1-AP



Basic Information

Catalog Number:

20218-1-AP

Size:

330 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003062

GeneID (NCBI):

6586

UNIPROT ID:

O75094

Full Name:

slit homolog 3 (Drosophila)

Calculated MW:

168 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

SLIT3, also named as KIAA0814, MEGF5 and SLIL2, may act as molecular guidance cue in cellular migration. It may be mediated by interaction with roundabout homolog receptors. The antibody is specific to SLIT3.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

