

For Research Use Only

NTN3-Specific Polyclonal antibody

Catalog Number: 20236-1-AP



Basic Information

Catalog Number:

20236-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006181

GeneID (NCBI):

4917

UNIPROT ID:

O00634

Full Name:

netrin 3

Calculated MW:

61 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

NTN3 also named as NTN2L. Netrins control guidance of CNS commissural axons and peripheral motor axons. This antibody is specific to NTN3.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

