

For Research Use Only

Mouse Trappc9, NIBP Polyclonal antibody, PBS Only

Catalog Number: 20287-1-PBS



Basic Information

Catalog Number:

20287-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14017

GenBank Accession Number:

BC070463

GeneID (NCBI):

76510

UNIPROT ID:

Q3U0M1

Full Name:

trafficking protein particle complex 9

Calculated MW:

129/250/128 kDa

Observed MW:

130-139 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20287-1-AP (Trappc9,NIBP antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 20287-1-PBS in a different storage buffer formulation.