

For Research Use Only

LEPROT Polyclonal antibody

Catalog Number: 20304-1-AP



Basic Information

Catalog Number:

20304-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14126

GenBank Accession Number:

BC011027

GeneID (NCBI):

54741

UNIPROT ID:

O15243

Full Name:

leptin receptor overlapping transcript

Calculated MW:

131 aa, 14 kDa

Observed MW:

16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, pig

Positive Controls:

WB : pig heart tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig heart tissue was subjected to SDS PAGE followed by western blot with 20304-1-AP (LEPROT antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.